



Non-Contact

FOREHEAD INFRARED THERMOMETER



0.5s
fast response



Body or surface
temperature mode



Clear LCD display

Medical Grade

DT-8806S

BYD CARES 



Display resolution
0.1°C/0.1°F



Memorization of the last
32 measurements



Backlight LCD display



Compact design

FEATURES

- ⊕ Precise non-contact measurements
- ⊕ Selectable °C or °F
- ⊕ Adjustable for body or surface temperature
- ⊕ Programmable alarm
- ⊕ Automatic data hold & Auto power off
- ⊕ Stores the last measurements up to 32 data points
- ⊕ Resolution 0.1°C(0.1°F)
- ⊕ Backlight LCD display

SPECIFICATIONS

Temperature range (body)	32.0°C to 42.5°C/89.6 to 108.5°F	
Temperature range (surface)	0°C to 60°C/32 to 140°F	
Display resolution	0.1°C/0.1°F	
Accuracy	32 to 35.9°C/89.6 to 96.6°F	±0.3°C/±0.5°F
ASTM	36 to 39°C/96.8 to 102.2°F	±0.2°C/±0.4°F
E1965-1998(2003)	39 to 42.5°C/102.2 to 108.5°F	±0.3°C/±0.5°F
Measuring distance	1-10cm/0.39-3.93in	
Response time	0.5s	

PRODUCT INFORMATION

Size	128x74x36mm / 5x2.9x1.4in
Weight	Net 130g
Packing/Parts	Manual , Warranty , Battery (2x1.5V “AAA”), Giftbox



BYD Signal & Communication Company Limited

Address : No.3009, BYD Road, Jiangling Community, Maluan Street, Pingshan District,
Shenzhen, Guangdong Province, P.R. China

Tel: +86 755 8988 8888 ext. 61673